



**TACMED**  
**SOLUTIONS**  
TACMEDSOLUTIONS.COM



 **TACMED**  
**MEDICAL**

# CONTROL WRAP

DESIGNED TO MAKE LIFE A LITTLE EASIER IN THE STRESSFUL SITUATION OF  
HEMORRHAGE CONTROL.

## PRODUCT SPECIFICATIONS

PRODUCT #		4" - CW4	6" - CW6
DIMENSIONS	4" Package	2"W x 4"L x 2"D   5.1cm W x 10.2cm L x 5.1cm D	
	6" Package	2"W x 6"L x 2"D   5.1cm W x 15.2cm L x 5.1cm D	
	4" Unrolled Length	72"L   182.9cm L	
	6" Unrolled Length	72"L   182.9cm L	
WEIGHT		6" Version - 2.2oz / 62g	4" Version - 1.6oz / 45g
NSN		6510-01-619-0260	



### INTENDED USE

To assist with device securement and hemorrhage control at the point of injury.



# CONTROL WRAP



Just as in our OLAES® and BLAST® Bandages, the **Control Wrap** includes a series of brakes that prevent accidental unrolling during application. The brakes also provide grip surfaces during the wrapping process and assist in maintaining the desired pressure and position of the bandage. This is extremely helpful when changing angles during the wrapping process or if bandaging awkward injuries.



## PRODUCT FEATURES

- 1 HOOK & LOOP CONTROL STRIPS**  
Built into the elastic wrap to prevent accidental unrolling during application.
- 2 SECURING CLIP**  
Locks the wrap into place to free up the user's hands to treat other injuries or casualties.



## ADDITIONAL FEATURES

### UNROLLING PREVENTION

By preventing accidental unrolling, you minimize the risks of contaminating a clean bandage, wasting a bandage due to contamination, losing time while managing loose ends, and losing pressure after applying the bandage due to shifting material.

### MATERIALS LIST

ABS, Elastic (Latex-Free), Poly 46 Thread.